

# PUR-VIEW™

**Disposable Flexible Bronchoscope/Videoscope**



**PUR-VIEW™** redefines bronchoscopy and airway management with its lightweight, ergonomic design and versatile sizing for both pediatric and adult patients. Engineered for procedures like **BAL**, **PDT**, and other critical interventions, it offers exceptional flexibility and HD imaging for confident navigation through complex anatomies. Designed to enhance safety, minimize contamination risk, and maximize procedural efficiency.

**PUR-VIEW™** - where clarity meets control.

**Sharper Images, Safer Procedures**

**MED-VERSE™**

# Key Features & Benefits



## Exceptional Visibility

High-definition vision with a premium resolution camera for accurate diagnosis and procedures.



## Unmatched Flexibility

210° up/down angulation enables controlled navigation, even in complex anatomical regions.



## Broad Compatibility

Designed for critical bronchoscopy applications in both routine and complex procedures.



## Instant Readiness

Plug & Play functionality ensures the scope is ready for immediate use, streamlining workflow.



## Enhanced Patient Safety

Single-use, disposable construction eliminates the risk of cross-contamination, making it ideal for high-risk patients.




## Seamless Connectivity

Equipped with USB and HDMI ports for fast, efficient data transfer and real-time viewing.



## Ergonomic Design

Lightweight (just 300g) and crafted for comfort, even during extended use.



**MED-VERSE™**



# Why PUR-VIEW™ ?

Superior  
Suction

Uncompromising  
Image Quality

Plug & Go  
Convenience

Enhanced User  
Experience



Suction Connector

Suction Button

Working Channel

**MED-VERSE™**



# ORDERING INFORMATION

|                     | Pur-View™   | Pur-View™   | Pur-View™   | Pur-View™   | Pur-View™   | Pur-View™   |
|---------------------|-------------|-------------|-------------|-------------|-------------|-------------|
| Models              | Mega        | Regular     | Regular     | Slim        | Slim        | UltraSlim   |
| Product Code        | DBR – MM581 | DBR – MM501 | DBR – MM401 | DBR – MM321 | DBR – MM251 | DBR – MM221 |
| Outer Diameter(mm)  | 5.8         | 5           | 4           | 3.1         | 2.5         | 2.2         |
| Working Channel(mm) | 3.4         | 2.8         | 2           | 1.2         | 1.2         | 0           |
| Max (mm)            | 6.2         | 5.5         | 4.5         | 3.2         | 2.5         | 2.3         |
| ETT (mm)            | 7           | 6           | 5           | 3.5         | 3           | 2.5         |
| Bending angle (°)   | 210/210     | 210/210     | 210/210     | 210/210     | 210/210     | 210/210     |
| FOV (°)             | 120         | 120         | 120         | 120         | 120         | 120         |
| DOV (mm)            | 3-50        | 3-50        | 3-50        | 3-50        | 3-50        | 3-50        |
| Working length (mm) | 600         | 600         | 600         | 600         | 600         | 600         |
| CMOS (pixel)        | 1280*720    | 1280*720    | 1280*720    | 1280*720    | 1280*720    | 1280*720    |
| Weight              | 300 gm      | 300 gm      | 300 gm      | 300 gm      | 300 gm      | 300 gm      |

|       |                   |               |
|-------|-------------------|---------------|
| MM-16 | Purview Processor | 16inch Screen |
|-------|-------------------|---------------|



**FOR MORE INFORMATION  
CONTACT US AT:**



0124-4712471



info@medorah.com



www.medorah.com

**MED-VERSE™**